APTOGLOBINA

Turbidimetria

REF. 6709 40+10 ml

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

USO PREVISTO

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BIBLIOGRAFIA

1. Clinical Guide to Laboratory Tests, Edited by NW Tietz W.B. Saunders Co., Philadelphia, 483 (1983).
2. Dati F. et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520
3. Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, (1987).
4. Kluthe R et al. Nature 205: 93-94 (1965).
5. Tarukosky PH et al. Scand J Clin Lab Invest 18:80-86 (1966).
6. Young DS. Effects of drugs on Clinical Laboratory Tests, 4th ed. AAC Press (1995).
7. Friedman and Young. Effects of disease on Clinical Laboratory tests, 3th ed. AAC Press (1997).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com - web site: www.giessediagnostics.com

670907

Ed. 2014/12 rev. 01

HAPTOGLOBIN

Turbidimetry

REF. 6709 40+10 ml

DNV CERTIFIED COMPANY

UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of Haptoglobin in serum and plasma.

PRINCIPLE

Anti-haptoglobin antibodies when mixed with samples containing haptoglobin, form insoluble complexes. These complexes cause an absorbance change, dependent upon the haptoglobin concentration of the patient sample, that can be quantified by comparison from a calibrator of known haptoglobin concentration.

SAMPLE

Fresh serum, plasma with EDTA or heparin.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolyzed or lipemic samples.

The samples are stable 7 days at 2-8°C or 3 months at -20°C.

KIT COMPONENTS

Reagent (A) Diluent Volume = 40 ml	Tris buffer 20 mmol/l, PEG 8000, pH 8.3 Sodium azide 0.95 g/l
Reagent (B) Antiserum Volume = 10 ml	Goat serum, anti-human haptoglobin, pH 7.5 Sodium azide 0.95 g/l

Optional: General Proteins Calibrator – REF. 7779

The Calibrator is not included in the kit.

Calibrator Volume = 2 ml	General Proteins Calibrator REF. 7779
-----------------------------	--

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C. Once opened, the reagents are stable at least 3 weeks at 2-8°C protected from light and in the absence of contamination. keep bottles closed when not in use.

REAGENT PREPARATION

The Reagents are ready to use.

Calibration Curve: Prepare the following dilutions of the Calibrator using NaCl 9 g/l as diluent.

Multiply the concentration of the haptoglobin calibrator by the corresponding factor stated in the below table to obtain the haptoglobin concentration of the each dilution.

Cal Dilution	1	2	3	4	5	6
Calibrator (µl)	--	10	25	50	75	100
NaCl 9 g/l (µl)	100	90	75	50	25	--
Factor	0	0.1	0.25	0.5	0.75	1.0

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength: 340 nm
Lightpath: 1 cm
Temperature: 37°C

Adjust the instruments to zero with distilled water

	sample	calibrator
pipette:		
Reagent (A)	800 µl	800 µl
sample	10 µl	
calibrator		10 µl

Mix and read the absorbance A1 after the sample addition.

Immediately add:

	200 µl	200 µl
Reagent (B)		

Mix and read the absorbance A2 of calibrators and sample exactly 2 minutes after the Reagent (B) addition.

RESULTS CALCULATIONS

Calculate the absorbance difference (A2 – A1) of each point of the calibration curve and plot the values obtained against the haptoglobin concentration of each calibrator dilution. Haptoglobin concentration in the sample is calculated by interpolation of its (A2 – A1) in the calibration curve.

EXPECTED VALUES

30 – 200 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control serum:

General Protein Control (REF. 7767).

PROCEDURE

Sensitivity: the sensitivity of the method is: 1.3 mg/dl. Values less than 1.3 mg/dl give non-reproducible results.

Prozone Effect No prozone effect up to 1200 mg/dl.

Linearity: the method is linear up to 300 mg/dl. For higher values, dilute the sample 1:5 and multiply the result by 5.

Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (mg/dl)	39.25	97.35	191.5
CV %	1.5	0.9	1.2

Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (mg/dl)	39.25	97.35	191.5
CV %	6.7	2.3	1.2

Interferences: bilirubin does not interfere up to 50 mg/dl. Hemoglobin up to 50 g/l and Rheumatoid factors up to 950 U/ml do not interfere. Lipemia in concentrations ≥ 6 g/l interfere

Correlation against a reference method: $Y = 0.88x + 4.8$ $r = 0.98$

REFERENCES

- Clinical Guide to Laboratory Tests, Edited by NW Tietz W.B. Saunders Co., Philadelphia, 483 (1983).
- Dati F. et al. Eur J Clin Chem Clin Biochem 1996; 34:517-520
- Pesce AJ and Kaplan, LA. Methods in Clinical Chemistry. The CV Mosby Company, St. Louis MO, (1987).
- Kluthe R et al. Nature 205: 93-94 (1965).
- Tarukosky PH et al. Scand J Clin Lab Invest 18:80-86 (1966).
- Young DS. Effects of drugs on Clinical Laboratory Tests, 4th ed. AAC Press (1995).
- Friedman and Young. Effects of disease on Clinical Laboratory tests, 3th ed. AAC Press (1997).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

670907

Ed. 2014/12 rev. 01