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Turbidimetria al lattice

REF. 6738 – 30+10

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

USO PREVISTO

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DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UNI EN ISO 13485:2012

INTENDED USE

Quantitative determination of glycated hemoglobin (HbA_{1c}) in blood.

PRINCIPLE

This method utilizes the interaction of antigen and antibody to directly determine the HbA_{1c} in whole blood. Total hemoglobin and HbA_{1c} have the same unspecific absorption rate to latex particles. When mouse antihuman HbA_{1c} monoclonal antibody is added (Reagent B) latex-antibody complex is formed.

Agglutination is formed when goat anti-mouse IgG polyclonal antibody interacts with the monoclonal antibody. The amount of agglutination is proportional to the amount of HbA_{1c} absorbed on the surface of latex particles. The amount of agglutination is measured as absorbance.

The HbA_{1c} value is obtained from a calibration curve.

SAMPLE

Whole blood with EDTA.

To determine HbA_{1c}, a hemolysate must be prepared for each sample:

- Dispense 1 ml Hemolysis Reagent (Reagent C) into tubes.
- Place 20 µl of well mixed whole blood (sample, calibrator or control). Mix.
- Allow to stand for 5 minutes or until complete lysis is evident.

KIT COMPONENTS

Reagent (A) Volume = 30 ml	Latex 0.13%, Buffer, Stabilizer
Reagent (B) Volume = 10 ml	Mouse anti-human HbA _{1c} monoclonal antibody 0.05 mg/ml, goat anti-mouse IgG polyclonal antibody 0.08 mg/ml. Stabilizers.
Reagent (C) Hemolysis reagent Volume = 2x100 ml	Water and stabilizers
Optional	Calibrator HbA _{1c} – REF. 6739 Control HbA _{1c} – REF. 6744

The calibrator and controls are not included in the kit.

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C. Do not use the reagents past their expiration date.

Once opened, the reagents are stable at least 2 months at 2-8°C protected from light and in the absence of contamination.

Keep bottles closed when not in use.

REAGENT PREPARATION

The reagents are ready to use.

Mix gently avoiding foam formation.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength: 630 - 660 nm

Lightpath : 1 cm

Temperature: 37°C

Adjust the instrument to zero with distilled water

pipette:

Reagent (A) 360 µl

Calibrators or sample 10 µl

Mix, incubate 5 minutes. Add:

Reagent (B) 120 µl

Mix, read the absorbance (A) after 5 minutes of the Reagent (B) addition.

RESULTS CALCULATION

HbA_{1c} concentration (%):

Plot (A) obtained against the HbA_{1c} concentration of each calibrator (1 to 4 Level).

HbA_{1c} percentage in the sample is calculated by interpolation of its (A) in the calibration curve.

EXPECTED VALUES

Non diabetics: < 6 %
Controlled diabetics < 7 %

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

For this purpose we recommend the use of control: **HbA_{1c} Control (REF 6744)**.

PERFORMANCE

Sensitivity: 1.0 % di HbA_{1c}.

Linearity: Range: 2.0 % - 16 %

Precision intra-assay:

	Level 1	Level 2
Mean	5.48	10.28
DS	0.078	0.176
CV %	1.43	1.72

Precision inter-assay:

	Level 1	Level 2
Mean	5.48	10.28
DS	0.152	0.275
CV %	2.77	2.68

Interferences:

1. Bilirubin up to 50 mg/dl, triglycerides up to 2000 mg/dl and ascorbic acid up to 50 mg/dl, do not interfere.
2. Results may be inconsistent in patients who have the following conditions: opiate addiction, lead-poisoning, alcoholism, ingest large doses of aspirin.

Correlation against a reference method: $Y = 1.050x - 0.481$ $r = 0.988$

REFERENCES

1. Trivelli, L.A., Ranney, H.M., and Lai, H.T., New Eng. J. Med. 284,353 (1971).
2. Gonen, B., and Rubenstein, A.H., Diabetologia 15, 1 (1978).
3. Gabbay, K.H., Hasty, K., Breslow, J.I., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J. Clin. Endocrinol. Metab. 44, 859 (1977).
4. Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978).
5. Tietz, N.W., Textbook of Clinical Chemistry, Philadelphia, W.B. Saunders Company, p. 794-795 (1999).
6. Ceriello, A. et al, Diabetologia 22, p. 379 (1982).
7. Little, R.R., et al, Clin Chem. 32, pp 358-360 (1986).
8. Fluckiger, R., et al, New Eng. J. Med. 304 pp 823-827 (1981).
9. Nathan, D.M., et al. Clin. Chem. 29, pp 93-97 (1989).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

673807

Ed. 2014/12 rev. 02