FERRITINA

Turbidimetria al Lattice

REF. 6719 40+10 ml

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

USO PREVISTO

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DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of Ferritin in serum.

PRINCIPLE

Latex particles coated with specific anti-human ferritin are agglutinated when they react with samples containing ferritin. The latex particles agglutination is proportional to the concentration of ferritin in the sample and can be measured by turbidimetry.

SAMPLE

Fresh serum. Stable 7 days at 2-8°C or 3 months at -20°C.

The samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolized or lipemic samples.

KIT COMPONENTS

Reagent (A) Diluent Volume = 40 ml	Glycine buffer 20 mmol/l, pH 8.5 Sodium azide 0.95 g/l
Reagent (B) Latex Volume = 10 ml	Latex particles coated with polyclonal anti-human ferritin antibodies, pH 8.2 . Sodium azide 0.95 g/l
Calibrator Volume = 3 ml	Human Ferritin The concentration is stated on the vial

Optional: Ferritin Control – REF. 7795

The control is not included in the kit.

Control Volume = 2 ml	Ferritin Control REF. 7795
--------------------------	-----------------------------------

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C . Once opened, the reagents are stable at least 3 weeks at 2-8°C protected from light and in the absence of contamination.

keep bottles closed when not in use.

Do not use the reagents after the expiration date.

REAGENT PREPARATION

The Reagents and the Calibrator are ready to use.

Calibration Curve: Prepare dilutions of the Calibrator using NaCl 9 g/l as diluent. Multiply the concentration of the Calibrator by the corresponding factor indicated in the table below to obtain the ferritin concentration of each point of the curve.

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PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Components from human origin have been tested and found to be negative for the presence of HbsAg, HCV and antibody to HIV 1/2. However handle cautiously as potentially infectious.

Dispose of waste according to local laws.

PROCEDURE

Wavelength: 650 nm (630)
Lightpath: 1 cm
Temperature: 37°C

Adjust the instrument to zero with distilled water

pipette:	sample	calibrator
Reagent (A)	800 µl	800 µl
Reagent (B)	200 µl	200 µl
sample	100 µl	
calibrator		100 µl

Mix and read the absorbance immediately (A1) and after 8 minutes (A2) of the latex addition.

RESULTS CALCULATION

Calculate the absorbance difference (A2 – A1) of each point of the calibration curve and plot the values obtained against the ferritin concentration of each calibrator dilution. Ferritin concentration in the sample is calculated by interpolation of its (A2 – A1) in the calibration curve.

EXPECTED VALUES

Men: 30 – 220 µg/l

Women: 20 – 110 µg/l

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control serum:

Ferritin Control (REF. 7795).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 10.0 µg/l. Values under 10.0 µg/l give non-reproducible results.

Prozone effect: No prozone effect up to 4000 µg/l .

Linearity: the method is linear up to 300 µg/l. For higher values, dilute the sample 1:5 and multiply the result by 5.

Precision intra-assay:

	Level 1	Level 2
Mean (µg/l)	65	178
CV %	3.56	1.87

Precision inter-assay:

	Level 1	Level 2
Mean (µg/l)	65	178
CV %	5.16	2.9

Accuracy: Results obtained with this reagents did not show systematic differences when compared with commercial reagents of similar characteristics.

Interferences: bilirubin does not interfere up to 20 mg/dl. Hemoglobin up to 10 g/l and rheumatoid factor up to 600 U/ml do not interfere. Lipemia interferes.

REFERENCES

1. Newman DJ, Henneberry H, Price CP. Ann Clin Biochem 29: 122-42 (1992).
2. Bernard A, Lauwerys R. J Immunol Methods. 71:141-147 (1984).
3. Wiedermann, Jonetz-Mentzel. Eur J Clin Chem Clin Biochem 31:453-457 (1993).
4. Worwood M. Blood Reviews. 4:259-269 (1990).
5. Lipschitz D, Cook JD, Finch CA. The New England J of Med 290: 1213-1216 (1974).
6. Mazza J, Barr RM, McDonald JWD. Can Med Assoc J 119:884-886 (1978).
7. Tietz Textbook of Clinical Chemistry, 3rd Ed. Burtis CA, Ashwood ER. WB Saunders Co., (1999).
8. Young DS. Effects of drugs on clinical laboratory tests. 3rd Ed. AACCC Press. (1997).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

671907

Ed. 2014/02 rev. 02