MICROALBUMINA

Turbidimetria

REF. 6740 40+10 ml

AZIENDA CERTIFICATA DNV

UNI EN ISO 9001:2008

UNI EN ISO 13485:2012

USO PREVISTO

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BIBLIOGRAFIA

1. Feldt-Rasmussen B et al. J Diab Comp 1994; 8:137-145
2. Panuyiotou B N. Journal International Medical Research 1994; 22:181-201
3. Bar J et al. Diabetic Medicine 1995; 12:649-656.
4. Gilbert R E et al. Diabetic Medicine 1994; 11:636-645.
5. Medcalf E A et al. Clin Chem 1990; 36/3: 446-449
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCPress (1995).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giesseiagnostics.com - web site: www.giesseiagnostics.com

674007

Ed. 2014/12 rev. 01

MICROALBUMIN

Turbidimetry

REF. 6740 40+10 ml

DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of Microalbumin in urine.

PRINCIPLE

Latex particles coated with specific antibodies anti-human albumin are agglutinated when mixed with samples containing microalbumin. The agglutination causes an absorbance change, dependent upon the microalbumin contents of sample that can be quantified by comparison from a calibrator of known microalbumin concentration.

SAMPLE

Urine 24 hours or first morning urine specimen.

Adjust the pH at 7.0 with NaOH/HCl 1 mol/l.

Stable: 7 days at 2-8°C with addition of sodium azide 1 g/l.

Centrifuge urine samples before testing.

KIT COMPONENTS

Reagent (A) Diluent Volume = 40 ml	Glycine buffer 100 mmol/l – pH 10.0 Preservative
Reagent (B) Latex Volume = 10 ml	Latex particles coated goat IgG anti-human albumin, pH 8.2. Preservative
Calibrator Volume = 1 ml	Liquid Calibrator. Concentration is stated on the label

Optional: Microalbumin Control – REF. 6750

The Control is not included in the kit.

Control Volume = 2 ml	Microalbumin Control REF. 6750
--------------------------	---------------------------------------

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C. Once opened, the reagents are stable at least 3 weeks at 2-8°C protected from light and in the absence of contamination. keep bottles closed when not in use.

REAGENT PREPARATION

Mix gently the latex before use.

Working solution: prepare the necessary amount as follows:

1 ml Latex (Reagent B) + 4 ml Diluent (Reagent A)

The working solution (A+B) is stable one day at 2-8°C.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Components from human origin have been tested and found to be negative for the presence of HbsAg, HCV and antibody to HIV 1/2. However handle cautiously as potentially infectious.

Dispose of waste according to local laws.

PROCEDURE

Wavelength: 540 nm

Lightpath: 1 cm

Temperature: 37°C

Adjust the instrument to zero with distilled water

pipette:	sample	calibrator
Working reagent (A+B)	1000 µl	1000 µl
sample	7 µl	
calibrator		7 µl

Mix and read the absorbance immediately (A1) and after 2 minutes (A2) of the sample addition.

RESULTS CALCULATION

$\mu\text{ALB (mg/L)} = (A2 - A1)_{\text{Sample}} / (A2 - A1)_{\text{Calibrator}} \times \text{Calibrator Value}$

EXPECTED VALUES

Up to **30 mg/24h** urine specimen and **20 mg/l** in a first morning urine specimen.

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control serum:

Microalbumin Control (REF. 6750).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 2 mg/l. Values less than 2 mg/l give non-reproducible results.

Prozone effect: No prozone effect was detected upon 1000 mg/l.

Linearity: the method is linear up to 150 mg/l. For higher values, dilute 1:5 with saline the samples and multiply the result by 5.

Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (mg/l)	10.36	16.95	57.33
CV %	1.9	1.4	1.1

Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (mg/l)	10.36	16.95	57.33
CV %	4.1	2.7	2.3

Interferences: bilirubin up to 10 mg/dl does not interfere. Hemoglobin does not interfere up to 10 g/l. Glucose does not interfere up to 2 g/l. Creatinine does not interfere up to 3 g/l. Urea in a concentration ≥ 1 g/l interfere.

Correlation against a reference method: $Y = 0.424x + 10.55$ $r = 0.99$

REFERENCES

1. Feldt-Rasmussen B et al. J Diab Comp 1994; 8:137-145
2. Panuyiotou B N. Journal International Medical Research 1994; 22:181-201
3. Bar J et al. Diabetic Medicine 1995; 12:649-656.
4. Gilbert R E et al. Diabetic Medicine 1994; 11:636-645.
5. Medcalf E A et al. Clin Chem 1990; 36/3: 446-449
6. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCPress (1995).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

674007

Ed. 2014/12 rev. 01