PCR
Turbidimetria con antisiero

REF. 6753 50+ 5 ml
REF. 6754 100+10 ml

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

USO PREVISTO

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DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of C-Reactive Protein (CRP) in human serum by turbidimetric immunoassay.

PRINCIPLE

Anti-human CRP antibodies when mixed with samples containing CRP, form a reaction antigen-antibody. This determines a turbidity measurable at 340 nm.

SAMPLE

Use fresh serum. The serum may be stored at 2-8°C for 48 hours. If stored for a longer period, the sample should be frozen.

KIT COMPONENTS

Reagent (A) - CRP Buffer Volume = 50/100 ml	Goods buffer, Sodium azide 0.95 g/l
Reagent (B) - CRP Antiserum Volume = 5/10 ml	Polyclonal Goat serum anti-human CRP, Sodium azide 0.95 g/l

Optional: CRP Calibrator – REF. 7774

The Calibrator is not included in the kit.

Calibrator Volume = 1.0 ml	CRP Calibrator REF. 7774 - 1x1 ml
-------------------------------	-----------------------------------

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C. Once opened, the reagents are stable at least 2 months at 2-8°C protected from light and in the absence of contamination. Do not use over expiry date. Keep bottles closed when not in use.

REAGENT PREPARATION

Reagents are ready to use.

Calibration curve: Prepare serial 1:2 dilutions of the Calibrator using NaCl 9 g/l as diluent: 1 ; 1:2 ; 1:4 ; 1:8 ; 1:16.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength:	340 nm
Lightpath:	1 cm
Temperature:	37°C
Method:	End point

Sample, Control or Calibrator	20 µl
Add Reagent (A)	300 µl
Mix, incubate for 1 minute, read absorbance, then add:	
Reagent (B)	30 µl
Mix, incubate and read absorbance after 4 minutes.	

Validated on BT2000 plus

EXPECTED VALUES

PCR: < 8 mg/l

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

PERFORMANCE

Measuring Range: 0 - 20 mg/dl.

Detection Limit: 0.1 mg/dl

Sensitivity: 0.018 ABS units / concentration unit.

Precision %CV:

	Level 1	Level 2	Level 3
Intra-assay:	1.13	1.65	1.17
Inter-assay:	2.31	0.98	1.01

Interferences: No interferences for: Hemoglobin (250 mg/dl), Sodium citrate (1000 mg/dl), Heparin (50 mg/dl), Bilirubin (20 mg/dl), EDTA (5 mg/dl), Triglycerides (2500 mg/dl).

Correlation against a reference method: $Y = 0.9905x - 0.1234$ $r = 0.99$

REFERENCES

1. Manack, J.R. and Richards, CB., J. Immunol. 20,1019 (1971).
2. Ritchie, RF., J.Lab.Clin. Med. 70, 512 (1967).
3. Pepys MB. et al., Ann. NY Acad. Sci, 389, 459 (1982).