# ALFA AMILASI SL

Metodo cinetico (CNP3)

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REF. 5501 6x 10 ml  
REF. 5502 3x 50 ml  
REF. 5503 3x 10 ml

AZIENDA CERTIFICATA DNV

UNI EN ISO 9001:2008

UN EN ISO 13485:2012

## USO PREVISTO

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Ed. 2014/12 rev. 01

# ALPHA AMYLASE SL

Kinetic Method (CNP3)  
Liquid Reagent ready to use

REF. 5501 6x 10 ml  
REF. 5502 3x 50 ml  
REF. 5503 3x 10 ml

DNV CERTIFIED COMPANY

UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## INTENDED USE

Quantitative determination of  $\alpha$ -Amylase in serum and plasma.

## PRINCIPLE

The  $\alpha$ -Amylase hydrolyzes 2-chloro-4-nitrophenyl- $\alpha$ -D-maltotriose (CNP3) into 2-chloro-4-nitrophenyl- $\alpha$ -D-maltoside (CNP2), maltotriose (G3), glucose and 2-chloronitrophenol. The absorbance change in unit time measured at 405 nm is proportional to the enzyme activity in the sample.

## SAMPLE

Serum, heparinized plasma, urine.

Do not use other anticoagulants like EDTA, citrate and oxalate as they inhibit the enzyme. Avoid hemolyzed samples.

Dilute urine 1:3 with saline solution.

Amylase activity is stable one week at 15-25°C and some months at 2-8°C.

## KIT COMPONENTS

|                                      |                  |            |
|--------------------------------------|------------------|------------|
| Reagent (A) AMY<br>Volume = 10/50 ml | Good buffer pH 6 | 100 mmol/l |
|                                      | CNP3             | 4.1 mmol/l |
|                                      | Sodium chloride  | 10 mmol/l  |
|                                      | Calcium acetate  | 1 mmol/l   |
|                                      | Sodium azide     | 15 mmol/l  |

The reagent is stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not freeze. Once opened reagent is stable for 2 months at 2-8°C if contamination is avoided.

Keep bottles closed when not in use.

## REAGENT PREPARATION

Liquid Reagent, bring to room temperature (15-25°C) before use.

## PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

### Warning

**H302** Harmful if swallowed.

**EUH032** Contact with acid liberates very toxic gas.

**P270** Do not eat, drink or smoke when using this product..

**P301+P312** IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

**P330** Rinse mouth.

**P273** Avoid release to the environment.

**Contains:** sodium azide.

Reagents may contain some non-reactive and preservative components. Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

## PROCEDURE

|                 |                         |
|-----------------|-------------------------|
| Wavelength:     | 405 nm                  |
| Lightpath :     | 1 cm                    |
| Temperature:    | 37°C                    |
| Reading:        | against distilled water |
| Method:         | Increasing Kinetic      |
| Sample/Reagent: | 1/50                    |

pipette:

|             |              |
|-------------|--------------|
| Reagent (A) | 1000 $\mu$ l |
| sample      | 20 $\mu$ l   |

Mix, incubate at 37°C for 1 minute, read the initial absorbance against water. Perform 3 readings at 60 seconds intervals. Calculate the average value of the absorbance variations per minute ( $\Delta A$ /min).

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

## RESULTS CALCULATION

Perform calculation in units per litre, multiplying the  $\Delta A$ /min by the factor:  
activity in U/L:  $\Delta A$ /min  $\times$  3517

## EXPECTED VALUES

Serum/plasma: 25 – 110 U/L  
Urine: < 480 U/L

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

## PERFORMANCE

**Sensitivity:** the sensitivity of the method is: 1 U/L

**Linearity:** the method is linear up to 2350 U/L. For higher values, dilute the sample 1:2 and multiply the result by 2.

### Precision intra-assay:

|            | Level1 | Level 2 | Level 3 |
|------------|--------|---------|---------|
| Mean (U/l) | 42.90  | 184.2   | 497.5   |
| DS         | 0.64   | 2.25    | 5.64    |
| CV %       | 1.49   | 1.22    | 1.13    |

### Precision inter-assay:

|            | Level 1 | Level 2 | Level 3 |
|------------|---------|---------|---------|
| Mean (U/l) | 44.15   | 181.2   | 524.5   |
| DS         | 0.48    | 2.82    | 2.22    |
| CV %       | 1.09    | 1.56    | 0.42    |

**Interferences:** bilirubin does not interfere up to 30 mg/dl. Triglycerides up to 1500 mg/dl do not interfere. Hemoglobin up to 500 mg/dl does not interfere. Glucose up to 500 mg/dl does not interfere.

**Correlation against a reference method:**  $Y = 0.4582x + 1.089$   $r = 0.999$

## REFERENCES

1. Lorentz, K., Clin. Chem. Clin. Biochem. 17, 499 (1979).
2. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
3. Young, D.S., et al., Clin. Chem. 21:1D (1975).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) - web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

550107

Ed. 2014/12 rev. 01