# BILIRUBINA TOTALE

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Ed. 2015/05 rev. 02

# TOTAL BILIRUBIN

Colorimetric Method  
Liquid Reagents ready to use

REF. 0037/50 4x 50 ml  
REF. 0037/2 2x100 ml  
REF. 0037 4x100 ml

DNV CERTIFIED COMPANY  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## INTENDED USE

Quantitative determination of Total Bilirubin in the serum.

## PRINCIPLE

In an acid medium, total bilirubin reacts with diazotized sulphanilic acid to form a pink diazo compound (azobilirubin), whose intensity is proportional to the concentration of bilirubin present in the sample.

Direct bilirubin consists of glucuronic acid conjugated derivatives, is water soluble and it reacts directly. Total bilirubin is obtained by the presence of a solubilizing agent that cleaves the bond with albumin.

## SAMPLE

Not hemolyzed fresh serum.

Analyze samples within 2 hours after collection. Protect samples from light. Direct sunlight may cause up to 50% decrease in direct bilirubin within 1 hour. The samples are stable 12 hours at 2-8°C in the dark and one month at -20°C.

## KIT COMPONENTS

Reagent (A) BIL-T Volume = 50/100 ml	Sulphanilic Acid Hydrochloric acid Solubilizing agent	46 mmol/l 270 mmol/l 10 mmol/l
Reagent (B) BIL-T Volume = 25 ml	Sodium nitrite	143 mmol/l

The reagents are stable until the expiration date when stored at temperatures as indicated on the label and protected from light. Do not use over expiry date. Once opened reagents are stable for 2 months if contamination is avoided. Keep bottles closed when not in use.

## REAGENT PREPARATION

Liquid Reagents, ready to use.

For use as **Monoreagent**: combine the Reagent (B) and (A) in the ratio 1:51. Add 0.5 ml of Reagent (B) to 25 ml of Reagent (A).

The working solution (A+B) is stable 5 days in refrigerator (2-8 °C) in a dark bottle, protected from light.

## PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

Reagent (A):

**EUH208** Contains: SULFANILIC ACID  
May produce an allergic reaction.

Reagents may contain some non-reactive and preservative components. Use the normal precautions required in the laboratory. Dispose of waste according to local laws.

## PROCEDURE

Wavelength: 536 nm (520 - 550)

Lightpath: 1 cm

Temperature: 37°C

Reading: against blank reagent

### Use as monoreagent:

pipette:	blank	sample	calibrator
Reagent (A+B)	1500 µl	1500 µl	1500 µl
water	100 µl		
sample		100 µl	
calibrator			100 µl

### Use as bi-reagent:

pipette:	blank	sample	calibrator
Reagent (A)	1200 µl	1200 µl	1200 µl
water	100 µl		
sample		100 µl	
calibrator			100 µl
Reagent (B)	40 µl	40 µl	40 µl

Mix and read the absorbance of blank sample (Abx) in 20 seconds. Incubate for 5 minutes at 37°C and read the absorbance of the sample (Ax) and the calibrator (Ac).

Reaction volumes can be proportionally modified.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

## RESULTS CALCULATION WITH FACTOR (BI-REAGENT)

Total bilirubin mg/dl = (Ax - Abx) x 21 at 546 nm

## RESULTS CALCULATION WITH CALIBRATOR

Total bilirubin mg/dl = (Ax - Abx)/(Ac - Abc) x Calibratore Value

**Human proteic calibrator Ref. 6002/8 is advised (Optional).**

Factor to convert mg/dl in µmol/l = 17.1

## EXPECTED VALUES

adults: 0.2 - 1.2 mg/dl ( 3.4 - 20.5 µmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

## PERFORMANCE

**Sensitivity:** the sensitivity of the method is: 0.06 mg/dl.

**Linearity:** the method is linear up to 11 mg/dl (in the use as monoreagent) and up to 25 mg/dl (for the bi-reagent). For higher values, dilute with saline 1:2 the sample and multiply the result by 2.

### Precision intra-assay:

	Level 1	Level 2	Level 3
Mean (mg/dl)	0.58	1.03	6.00
DS	0.019	0.020	0.052
CV %	3.27	1.93	0.86

### Precision inter-assay:

	Level 1	Level 2	Level 3
Mean (mg/dl)	0.55	1.05	5.82
DS	0.011	0.029	0.104
CV %	1.98	2.75	1.78

**Interferences:** triglycerides up to 2000 mg/dl do not interfere.

Hemoglobin up to 150 mg/dl does not interfere.

**Correlation against a reference method:** Y = 0.975x + 0.042 r = 0.9994

## REFERENCES

1. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
2. J.A. Lott and B.T. Doumas, Clin. Chem. 641-647, 39 (1993).
3. Young D. S., et al, Clin. Chem. 21:1D (1975).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18,300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) - web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

003707

Ed. 2015/05 rev. 02