CLORURI

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Ed. 2014/12 rev. 01

CHLORIDE

Mercury Thiocyanate Colorimetric Method
Liquid Reagent ready to use

REF. 0020 2x100 ml
REF. 0049 4x100 ml

DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of Chlorides in serum, plasma and urine according to mercury thiocyanate reaction.

PRINCIPLE

The chloride ions present in the sample react with mercuric ions, releasing the same quantity of thiocyanate ions which, in the presence of ferric ions, forms a red colored complex. The color intensity is proportional to the chloride concentration present in the sample.

SAMPLE

Serum, heparinized plasma, urine 24/h diluted 1:2 with distilled water and acidified with 2-3 drops of HCl 23 %.

Chlorides in the sample are stable up to 7 days at 2-8°C and one month at -20°C.

KIT COMPONENTS

| | | |
|-----------------------------------|---------------------|------------------------|
| Reagent (A) Cl Volume = 100 ml | Mercury thiocyanate | 25 mmol/l |
| | Mercury nitrate | 0.2 mmol/l |
| | Iron nitrate | 25 mmol/l |
| | Nitric acid | 60 mmol/l |
| Standard Volume = 10 ml | Sodium Chloride | 100 mEq/l (100 mmol/l) |
| | Sodium azide | 14 mmol/l |

The reagents are stable until the expiration date indicated on the label if stored at 15-25°C and protected from light. Do not use over expiry date.

Once opened reagents are stable for 2 months if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagent ready to use.

PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

Reagent (A):

Danger

- H370** Causes damage to organs.
H332 Harmful if inhaled.
H312 Harmful in contact with skin.
H302 Harmful if swallowed.
P264 Wash thoroughly after handling.
P280 Wear protective gloves/ protective clothing / eye protection / face protection.
P301+P312 IF SWALLOWED: Call a POISON CENTER or doctor / physician if you feel unwell.
P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
P405 Store locked up.

Contains: Methanol.

Reagents may contain some non-reactive and preservative components. Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

| | |
|-----------------|-----------------------|
| Wavelength: | 480 nm (450 – 500) |
| Lightpath: | 1 cm |
| Temperature: | 25, 30, 37°C |
| Reading: | against blank reagent |
| Method: | Increasing End Point |
| Sample/Reagent: | 1/100 |

| | | | |
|-------------|---------|---------|----------|
| pipette: | blank | sample | standard |
| Reagent (A) | 2000 µl | 2000 µl | 2000 µl |
| water | 20 µl | | |
| sample | | 20 µl | |
| standard | | | 20 µl |

Mix, incubate at 25, 30, 37 °C for 5 minutes. Read the absorbance of the sample (Ax) and the standard (As) against the blank reagent.

The reaction volumes can be proportionally modified.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific application.

RESULTS CALCULATION

Serum, plasma:

$$\text{Chloride mEq/l} = \text{Ax/As} \times 100 \text{ (Standard Value)}$$

Urine 24/h:

$$\text{Chloride mEq/l} = \text{Ax/As} \times 100 \times 2 \text{ (dilution factor)} \times \text{Urine Vol. 24/h (L)}$$

EXPECTED VALUES

Serum/plasma: **95 – 110 mEq/L**

Urine: **160 – 250 mEq/24h**

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 1 mEq/L.

Linearity: the method is linear up to 140 mEq/L. For higher values, dilute the sample 1:2 with saline and multiply the result by 2.

Precision intra-assay:

| | Level 1 | Level 2 | Level 3 |
|--------------|---------|---------|---------|
| Mean (mEq/l) | 46.89 | 88.80 | 117.5 |
| DS | 0.25 | 0.29 | 0.53 |
| CV % | 0.53 | 0.33 | 0.45 |

Precision inter-assay:

| | Level 1 | Level 2 | Level 3 |
|--------------|---------|---------|---------|
| Mean (mEq/l) | 46.15 | 90.54 | 121.1 |
| DS | 0.22 | 0.55 | 0.74 |
| CV % | 0.48 | 0.61 | 0.61 |

Interferences: bilirubin does not interfere up to 20 mg/dl. Triglycerides do not interfere up to 300 mg/dl. For particularly icteric, hemolyzed and lipemic sample, prepare a blank sample by adding 2000 µl of saline solution to 20 µl of sample. Read at 480 nm the absorbance of the blank sample against a saline solution and then deduct it from the sample absorbance (Ax).

Correlation against a reference method: $Y = 0.9898x + 1.0063$ $r = 0.9871$

REFERENCES

1. Kaplan LA, Pesce AJ, Clinical Chemistry, Mosby (1989)
2. Shales O. and Shales S., J. Biol. Chem. 140, 879 (1941)
3. Young D.S. et al., Clin. Chem. 21:1D (1975)

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com - web site: www.giessediagnostics.com

04907

Ed. 2014/12 rev. 01