# NUMERO DI DIBUCAINA

Metodo cinetico (Butiriltilcolina)

REF. 0029 6x10 ml

AZIENDA CERTIFICATA DNV  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## USO PREVISTO

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002907  
Ed. 2014/11 rev. 01

# DIBUCAINE NUMBER

Kinetic Method (Butyrylthiocholine)

REF. 0029 6x10 ml

DNV CERTIFIED COMPANY  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## INTENDED USE

Quantitative determination of Dibucaine Number.

## PRINCIPLE

Cholinesterase catalyzes the hydrolysis of butyrylthiocholine, forming butyrate and thiocholine, which reduces the ferricyanide ions to ferrocyanide.

The decrease in absorbance is followed at 405 nm and it is proportional to cholinesterase activity in examined sample.

Cholinesterase activity is measured with and without presence of dibucaine as inhibitor. Dibucaine number is calculated on the basis of inhibition percentage.

## SAMPLE

Serum, plasma with heparin.

The activity in the sample is stable 15 days at 2-8°C.

## KIT COMPONENTS

Reagent Volume = 10 ml	Dibucaine	16 mmol/l
---------------------------	-----------	-----------

The reagent is stable until the expiration date indicated on the label if stored at 2-8°C and protected from light. Do not use over expiry date.

Once opened reagent is stable for 2 months at 2-8°C if contamination is avoided. Do not freeze.

Keep bottles closed when not in use.

## REAGENT PREPARATION

Liquid Reagent, bring to room temperature (15-25°C) before use.

The Reagent Dibucaine is used with the Kit CHOLINESTERASE SL (REF. 0016-4153).

Reagent Preparation with inhibitor:

REAGENT (A) CHOLINESTERASE SL	REAGENT (B) CHOLINESTERASE SL	DIBUCAINE
5 ml	1 ml	1,5 ml

## PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

## PROCEDURE

Wavelength:	405 nm
Lightpath :	1 cm
Temperature:	37°C
Reading:	against distilled water
Method:	Decreasing kinetic

pipette:	REACTION WITH INHIBITOR	REACTION WITHOUT INHIBITOR
Reagent (A+B) with Dibucaine	600 µl	---
Reagent(A+B) without Dibucaine	---	600 µl
Sample	10 µl	10 µl

Mix, incubate at 37°C for 1 minute, read the initial absorbance against water. Perform 3 readings at 30 seconds intervals. Calculate the average value of the absorbance variations per minute. ( $\Delta A/\text{min}$ ).

Reaction volumes can be proportionally varied.

## RESULTS CALCULATION

Perform calculation in Units per litre, multiplying the  $\Delta A/\text{min}$  by the factor as it is indicated:

CHE total activity U/L:	$\Delta A/\text{min} \times 65800 (*)$
CHE inhibited activity U/L:	$\Delta A/\text{min} \times 65800 (*)$

(\*) If you not use the Kit CHOLINESTERASE SL of Giesse Diagnostics (REF. 0016-4153), multiply by factor provided by the cholinesterase in use.

Dibucaine Number Calculation (DN):

$$DN = 100 - [ ( \text{CHE inhibited activity} / \text{CHE without inhibitor activity} ) \times 100]$$

## EXPECTED VALUES (DIBUCAINE)

Normal Subjects	70 – 90 %
Heterozygous	30 – 70 %
Homozygotes	0 – 20 %

## EXPECTED VALUES (CHOLINESTERASE)

Men:	5100 - 11200 U/L
Women:	4200 - 10800 U/L

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

## PERFORMANCE

**Sensitivity:** the sensitivity of the method is: 160 U/L

**Linearity:** the method is linear up to 25 KU/L. For higher values, dilute with saline 1:3 the sample and multiply the result by 3.

**Precision intra-assay:**

	Level 1	Level 2	Level 3
Mean (U/L)	3850	6735	13749
DS	57.04	155.22	166.17
CV %	1.46	2.30	1.21

**Precision inter-assay:**

	Level 1	Level 2	Level 3
Mean (U/L)	3850	6744	13758
DS	46.70	92.66	126.87
CV %	1.21	1.37	0.92

**Interferences:** bilirubin does not interfere up to 20 mg/dl. A moderate hemolysis does not interfere in the results.

**Correlation against a reference method:**  $Y = 0.3194x + 527$   $r = 0.9831$

## REFERENCES

1. Eur.J.Clin. Chem.Clin. Biochem. Vol. 30, 1992, 162 -170
2. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
3. Tietz Textbook of Clinical Chemistry, 2th Ed., Burtis-Ashwood (1994)

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I.V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) -Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) – web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

002907  
Ed. 2014/11 rev. 01