FRUTTOSAMINA SL

Metodo cinetico colorimetrico NBT

REF. 5531 4x10 + 1x10 ml

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

USO PREVISTO

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DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of fructosamine in serum.

PRINCIPLE

At alkaline pH values, nitro-blue tetrazolium (NBT) is reduced by the action of fructosamine to give rise to the production of formazane, whose rate is directly proportional to the serum fructosamine concentration.

SAMPLE

Fresh serum unhemolyzed.

Fructosamine in serum are stable 3 days at room temperature (15-25°C), 2 weeks at 2-8°C and 2 months at -20°C.

KIT COMPONENTS

Reagent (A) Volume = 4x10 ml	Buffer NBT	100 mmol/l 0.48 mmol/l
Reagent (B) Volume = 1x10 ml	Uricase	≥ 1000 U/l
Calibrator (C) Lyo for 1 ml	Fructosamine	Value on the label

The reagent is stable until the expiration date indicated on the label if stored at 2-8°C, in the dark. Do not use over expiry date.

Once opened reagent is stable for 2 months at 2-8°C if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagents, bring to room temperature before use.

Reconstitute the Calibrator (C) with 1 ml of distilled water, divide it in several tubes and frozen. Stability: 2 months at -20°C.

PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

PROCEDURE

Wavelength:	530 nm (530 – 550)
Lightpath :	1 cm
Temperature:	37°C
Reading:	against water
Method:	Increasing Kinetic

pipette:	Sample	Calibrator
Reagent (A)	1000 µl	1000 µl
Reagent (B)	250 µl	250 µl
sample	100 µl	
calibrator		100 µl

Mix, incubate at 37°C for 10 minutes and read against water the absorbance of the sample (Ax) and the calibrator (Ac). Repeat the reading after 4 minutes. Calculate the average value of the absorbance variations per minute. ($\Delta A/\text{min}$).

Reaction volumes can be proportionally varied.

This method describes the manual procedure to use the kit.

For automated procedure, ask for specific applications.

RESULTS CALCULATION

$$\text{Fructosamine } (\mu\text{mol/l}) = \Delta A_x / \Delta A_c \times \text{Calibrator Value}$$

EXPECTED VALUES

Fructosamine: **184 – 290 µmol/l**

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

PERFORMANCE

Sensitivity: the sensitivity of the method is 10 µmol/l.

Linearity: the method is linear up to 1000 µmol/l. For higher values, dilute the sample 1:2 and multiply the result by 2.

Precision intra-assay:

	Level 1	Level 2
Mean (µmol/l)	105	352
DS	1.22	1.20
CV %	2.60	2.24

Precision inter-assay:

	Level 1	Level 2
Mean (µmol/l)	102	348
DS	1.76	1.96
CV %	1.21	1.14

Interferences: bilirubin does not interfere up to 2 mg/dl. Hemoglobin does not interfere up to 100 mg/dl. Triglycerides do not interfere up to 500 mg/dl. Ascorbic acid does not interfere up to 3 mg/dl.

Condition of altered protein metabolism can affect the determination of fructosamine.

In hydremic conditions (for example pregnancy), it is advisable to relate fructosamines to proteins by the following formula:

$$\text{Rel. Fructosamine } (\mu\text{mol/l}) = \text{fructosamine} \times 7.2 \text{ (in } \mu\text{mol/l) / Tot Proteins (in g/dl)}$$

Correlation against a reference method: $Y = 1.02x - 0.06$ $r = 0.999$

REFERENCES

1. Johnson R.N., Metcalf P.a., Baker J.R., Clin.Chim.Acta 127, 87-95 (1983).
2. Schleicher E.D., Vogt B. W., Clin. Chem. 36, 136-139 (1990).
3. Young D. S., et al, Clin. Chem. 21:1D (1975).

Giessere Diagnostica srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

553107
Ed. 2018/03 rev. 06