OMOCISTEINA SL

Metodo enzimatico bireagente

REF. 5540	2x15 + 2x2 ml + Cal 2x1ml
REF. 5539	15 + 2 ml + Cal 2x1ml
REF. 5543	2x15 + 2x2 ml (senza Cal)
REF. 5538	15 + 2 ml (senza Cal)

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UNI EN ISO 13485:2012

SIGNIFICATO CLINICO

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HOMOCYSTEINE SL

Bireagent enzymatic method

REF. 5540	2x15 + 2x2 ml + Cal 2x1 ml
REF. 5539	15 + 2 ml + Cal 2x1 ml
REF. 5543	2x15 + 2x2 ml (without Cal)
REF. 5538	15 + 2 ml (without Cal)

DNV CERTIFIED COMPANY

UNI EN ISO 9001:2008

UN EN ISO 13485:2012

CLINICAL SIGNIFICANCE

Homocysteine (HCY) is a thiol-containing amino acid produced by the intracellular demethylation of methionine. HCY is exported into plasma where it circulates mostly in its oxidized form bound to plasma proteins as a protein-HCY mixed disulphide with albumin (SS-HCY proteic). Smaller amounts of reduced homocysteine and the disulphide homocysteine (HCY-SS-HCY) are present. Total homocysteine (tHCY) represents the sum of all HCY species found in plasma and serum (free plus protein bound). HCY is either metabolized to cysteine or to methionine. In the vitamin B6 dependent trans-sulphuration pathway HCY is irreversibly catabolized to cysteine. A major part of HCY is remethylated to methionine, mainly by the folate and cobalamin-dependent enzyme methionine synthase. HCY accumulates and is excreted into the blood when these reactions are impaired. Severely elevated concentrations of HCY are found in subjects with homocystinuria, a rare genetic disorder of the enzymes involved in the metabolism of HCY. Patients with homocystinuria exhibit mental retardation, early arteriosclerosis and arterial and venous thromboembolism. Other less severe genetic defects that lead to moderately elevated levels of HCY are also found. Patients with chronic renal disease experience an excess morbidity and mortality due to arteriosclerosis CVD. Elevated concentration of HCY is a frequently observed finding in the blood of these patients. Although such patients may lack some of the vitamins involved in the metabolism of HCY, the increased levels of HCY are mainly due to impaired removal of HCY from the blood by the kidney. Recently there have been evidence of abortions and baby disease because of an elevated concentration of HCY.

PRINCIPLE

Bound or dimerised HCY (oxidised form) is reduced to free HCY, which then reacts with serine catalysed by cystathionine beta-synthase (CBS) forming cystathionine. Cystathionine in turn is broken down by cystathionine beta-lyase (CBL) to form homocysteine, pyruvate and ammonia. Pyruvate is then converted by lactate dehydrogenase (LDH) to lactate with nicotinamide adenine dinucleotide (NADH) as coenzyme.

The rate of NADH conversion to NAD⁺ is directly proportional to the concentration of homocysteine (ΔA 340 nm).

Reduction: Dimerised homocysteine, mixed disulfide, and protein-bound forms of HCY in the sample are reduced to form free HCY by the use of tris [2-carboxyethyl] phosphine (TCEP).

Enzymatic Conversion: Free HCY is converted to cystathionine by the use of cystathionine beta-synthase and excess serine. The cystathionine is then broken down to homocysteine, pyruvate and ammonia. Pyruvate is converted to lactate via lactate dehydrogenase with NADH as coenzyme. The rate of NADH conversion to NAD⁺ (ΔA 340 nm) is directly proportional to the concentration of homocysteine.

SAMPLE

Serum (collected in serum or serum separator tubes) and plasma (collected in potassium EDTA or lithium heparin tubes) may be used for the measurement of homocysteine. However, it is not recommended to use individual patient results from serum, heparinized plasma and EDTA plasma interchangeably. Additionally matrix differences between serum and serum separator tubes and plasma tubes have been reported.

To minimize increases in homocysteine concentration from synthesis by red blood cells process specimens as follows:

- Place all specimens (serum and plasma) on ice after collection and prior to processing. Serum may clot more slowly and the volume may be reduced.
- All specimens may be kept on ice for up to 6 hours prior to separation by centrifugation.
- Separate red blood cells from serum or plasma by centrifugation and transfer to a sample cup or other clean container.

Note: 1. Specimens not placed on ice immediately may exhibit a 10-20% increase in homocysteine concentration.

2. If the assay will be performed within 2 weeks after collection, the specimen should be stored at 2-8°C. If the testing will be delayed more than 2 weeks, the specimen should be stored frozen at -20°C or colder.

Specimens have been shown to be stable at -20°C for 8 months.

3. It is the responsibility of the operator to verify the correct specimen type(s) (are) used with the Homocysteine LS reagent.

4. Inspect all samples (specimens, calibrators and controls) for bubbles. Remove bubbles prior to analysis.

Giessa Diagnostics srl

Via Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM)

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com - web site: www.giessediagnostics.com

5. Specimens containing particulate matter (fibrin, red blood cells, or other matter) and visibly lipemic specimens should not be used with the assay. Results from these specimens may be inaccurate.

6. Mix specimens **thoroughly** after thawing by low speed vortexing or by gentle inversion to ensure consistency in results. Avoid repeated freezing and thawing. Specimens showing particulate matter, erythrocytes, or turbidity should be centrifuged before testing.

7. On board instrument storage. Plasma samples can be stored for 3 hours on-board the instrument.

PACKAGE

R.(A). Colourless and odourless liquid 2x15mL (Ref. 5540) or 1x15 mL (Ref. 5539) in amber vial NADH (0,47 mM), LDH (38 KU/L), Serine (0,76 mM), Trizma Base 1-10 %, Trizma hydrochloride 1-10%, Sodium azide < 1%. Reductant (TCEP:2,9 mM) **Ready to use.**

R.(B). Pale yellow odourless liquid 2x2,0 mL (Ref. 5540) or 1x2.0 mL (Ref. 5539) in amber vial. Cycling Enzymes CBS (0,748 KU/L) and CBL (16,4 KU/L) Sodium azide < 1%. **Ready to use.**

CAL LEVEL I. Colourless and odourless liquid 1x1.0 mL in vial. Aqueous Homocysteine blank (0 µmol/L). **Ready to use.**

CAL LEVEL II. Colourless and odourless liquid 1x1.0 mL in vial. Aqueous Homocysteine solution (28 µmol/L). **Ready to use.**

Calibrators are prepared gravimetrically and are traceable to NIST SRM 1955, confirmed by a designated measurement procedure (HPLC). The assigned values are printed on the labels (0 µmol/L e 28 µmol/L).

EXPECTED VALUES

Reference Range: The reference range should be determined by each laboratory to confirm the characteristics of the population being tested. As a point of reference the following data may be used until the laboratory has analysed a sufficient number of specimens to determine its own reference range. The HCY concentration in plasma or serum of healthy individuals varies with age, gender, geographical area and genetic factors. Scientific literature reports reference values for adult male and females between **5 and 15 µmol/L**, men having higher values than women, and post menopausal woman having higher homocysteine values than pre-menopausal women. HCY values will normally increase with age, giving a reference range among an elderly population (> 60 years) of 5-20 µmol/L. In countries with folic acid fortification programmes, reduced levels of HCY may be observed.

Measurable Range: The measurable range of the homocysteine assay with kit 5540 (or 5539) is 1-46 µmol/L.

WASTE DISPOSAL

Apply the rules laid down in D.Leg.vo 22/97 and subsequent amendments (Special and Hazardous Waste with relevant CER code).

STABILITY AND STORAGE

Store Reagents and calibrators at 2-8°C. **DO NOT FREEZE.** Reagents and calibrators may be used until the expiry date on the labels if stored as directed. Do not mix different reagent kit lot number. Reagents and calibrators are liquid and ready to use.

PROCEDURE

Wavelength:	340 nm
Optical path:	1 cm
Temperature:	37 °C
Reading:	against blank reagent
Reaction:	Fixed time
Calibration:	Multipunti

	Blank	Sample	Calibrator
Reagente (A)	280 µl	280 µl	280 µl
Water	20 µl		
Sample		20 µl	
Calibrators (1-2)			20 µl
Mix and incubate for 4 minutes			
Reagent (B)	28 µl	28 µl	28 µl

Shake and put in cuvette and, after the incubation time of 1 minutes at 37 ° (T.O'), read the calorimeters and sample extinction at 0" and after 220". Calculate the $\Delta E/\text{min}$ at 340 nm.

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On OLYMPUS AU400 the calibration curve is stable for 30 days.

PERFORMANCE

A correlation study was performed with plasma specimens from apparently healthy adults. All specimens were analysed using the homocysteine LS kit according to the CLSI (formally NCCLS) document EP9-A2. All results are based on 95% confidence Interval. All homocysteine kits specimens ranged from 6.5 to 49.0 $\mu\text{mol/L}$. The data obtained gave the following statistical values:

Comparison Method Number of specimens 94 Slope of regression line 0.991 Y-Intercept 0.165

Correlation coefficient 1.0

Precision

A study was performed with guidance from the CLSI (formally NCCLS) Document EP5-A2.

Three HCY controls and three human plasma panels were assayed using two lots of reagents, in replicants of two, at two separate times per day for 20 days on one instrument (n=80). A calibration curve was generated at the start of the study and was used throughout.

Dilution Linearity

The dilution linearity of the specimens gives a % recovery range of 91–104% for all samples across the range of the assay (1–46 $\mu\text{mol/L}$) on the OLYMPUS AU400.

Samples > 46 $\mu\text{mol/L}$ exhibit mean recovery of 100% + 11% of the expected result when diluted into the assay range.

Limit of Detection

The limit of detection (LOD) of the homocysteine kit according to the CLSI (formally NCCLS) Document EP17-A 30 was found to be 0.33 $\mu\text{mol/L}$.

Analytical Specificity

The specificity of the homocysteine LS was assessed according to guidance in the CLSI Document EP7-A231 for the interfering substances listed in the table below:

Bilirubin 20 mg/dl $\leq +10$	Glutathione 1000 $\mu\text{mol/l} \leq +10$
Haemoglobine 500 mg/dl $\leq +10$	Methionina 800 $\mu\text{mol/l} \leq +10$
RBC 0.4 % $\leq +10$	Cysteina 200 $\mu\text{mol/l} \leq +10$
Triglycerides 500 mg/dl $\leq +10$	Pyruvate 1250 $\mu\text{mol/l} \leq +10$

None of these substances interfered significantly in the assay.

Samples with raised protein levels show > 10% difference compared to results obtained from normal samples and should be avoided.

METHOD LIMITATIONS

- Samples with a HCY value higher than 50 $\mu\text{mol/l}$ should be diluted
- The reagent should be clear. Discard if turbid or if the beginning absorption is less than 0.5 to 340 nm .
- The S-adenoyilhomocysteine (SAH) may cause a remarkable and positive interference.
- Patients who took the following drugs (methorexato, carbamazepine, phenytoin, nitrous oxide, anticonvulsive,6-azuridine triacetate), may have high level of Hcy caused by the interference between drugs and HCY metabolism.
- It has been suggested to add 3-deazaadenosine to inhibit the production of Hcy in the red cells. It is not possible the use of a sample containing 3 deazaadenosine as this inhibits one of the enzyme used for the production of the kit.

QUALITY CONTROL

It is recommended to execute the quality control at every kit utilization to verify that values are within the reference range indicated by the methodology. For this purpose the use of test serum REF.5541 is advised.

PRECAUTIONS AND WARNINGS

Each time infected agents, chemical reagents, human or animal reagents, blood and every kind of biological substance are handled with, it is recommended to follow the common hygienic precaution; ex: wearing gloves, etc.

Safety Precautions

1. Adhere strictly to the instructions in this leaflet, particularly for handling and storage conditions.
2. Reagent A and Reagent B contain sodium azide which can react with lead or copper plumbing to form highly explosive metal azides. On disposal, flush with large quantities of water .
3. Material safety data sheets for all hazardous components contained in this kit are available upon request.

LIMITATIONS

1. The linear range of the homocysteine kit when run as directed is 1-46 $\mu\text{mol/L}$. Specimens > 46 $\mu\text{mol/L}$ should be diluted 1 part specimen to 2 parts Cal 0 $\mu\text{mol/L}$ or 1 part specimen to 9 parts Cal 0 $\mu\text{mol/L}$ as appropriate.
2. The Reagents should be clear.
3. Cystathionine is measured with homocysteine, but in the general population the cystathionine level (0.065 to 0.3 $\mu\text{mol/L}$) has a negligible effect. In very rare cases, renal disease and patients with severe metabolic disturbances, cystathionine levels may rise dramatically and in severe cases cause greater than 20% interference.
4. Hydroxylamine, present in several iron reagents may carryover (reagent probe or reaction cuvette) and cause falsely low results. Routine rinsing procedures are not adequate to eliminate this problem in most cases. Possible solutions would include special washing protocols, changing to an iron assay that used ascorbic acid as reductant or running iron and homocysteine assays on separate instruments.
5. Carbamazepine, methotrexate, phenytoin, nitrous oxide, or 6-azauridine triacetate may affect the homocysteine concentration.
6. Samples with raised protein levels show > 10% difference compared to results obtained from normal samples and should be avoided.
7. Note: Specimens from patients who are on drug therapy involving S-adenosyl-methionine may show falsely elevated levels of homocysteine. Patients who are taking methotrexate, carbamazepine, phenytoin, nitrous oxide, anticonvulsants, or 6-azauridine triacetate, may have elevated levels of homocysteine due to their effect on the pathway.
8. Specimens containing particulate matter (fibrin, red blood cells, or other matter) and visibly lipemic specimens should not be used with the assay. Results from these specimens may be inaccurate.

REFERENCES

1. McCully KS. Vascular Pathology of Homocysteinemia: Implications for the Pathogenesis of Arteriosclerosis. *Am J Pathol* 1969;56:111-122.
2. Malinow MR. Plasma Homocyst(e)ine and Arterial Occlusive Diseases: A Mini-Review. *Clin Chem* 1995;41:173-176.
3. Ueland PM. Homocysteine Species as Components of Plasma Redox Thiol Status. *Clin Chem* 1995;41:340-342.
4. Finkelstein JD. Methionine Metabolism in Mammals. *J Nutr Biochem* 1990; 1:228-237.
5. Mudd SH, Levy HL, Skovby F. Disorders of Transsulfuration. In: Scriver CR, Beaudet AL, Sly WS, et al., eds *The Metabolic and Molecular Basis of Inherited Disease*. New York: McGraw- Hill, 1995; 1279-1327.
6. Clarke R, Daly L, Robinson K, et al. Hyperhomocysteinemia: An Independent Risk Factor for Vascular Disease. *N Engl J Med* 1991;324:1149-1155.