# LATTATO

Metodo enzimatico colorimetrico

REF. 6751 10x10 ml  
REF. 6752 5x10 ml

AZIENDA CERTIFICATA DNV  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## USO PREVISTO

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AZIENDA CERTIFICATA DNV  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## INTENDED USE

Quantitative determination of lactate in serum and plasma.

## PRINCIPLE

Lactate is oxidized by lactate oxidase (LO) to pyruvate and hydrogen peroxide, which under the influence of peroxidase (POD), 4-aminophenazone (4-AP) and 4-chlorophenol forms a red quinone compound. The intensity of the color formed is proportional to the lactate concentration in the sample.

## SAMPLE

Serum, plasma with heparin. Do not use hemolyzed samples.  
Serum or plasma must be separated of the blood cells within 15 minutes; the reason is that blood cells will metabolise glucose to lactic acid.

## KIT COMPONENTS

Reagent (A) Buffer Volume = 100 ml	PIPES pH 7.5 4-chlorophenol	50 mmol/l 4 mmol/l
Reagent (B) Enzymes	Lactate oxidase Peroxidase 4-aminophenazone	800 U/l 2000 U/l 0.4 mmol/l
Calibrator	Lactate aqueous primary Std	10 mg/dl

The reagents are stable until the expiration date indicated on the label if stored tightly closed at 2-8°C and protected from light.

Once opened, the calibrator is stable one month at 2-8°C, tightly closed, protected from light, in the absence of contamination.

**Signs of reagent deterioration:** Presence of particles or turbidity.

Blank absorbance (A) at 580 nm  $\geq 0.18$

## REAGENT PREPARATION

Dissolve the contents of one vial of Reagent (B) in 10 ml of Reagent (A).  
Mix gently.

The working solution (A+B) is stable one month at 2-8°C or one week at room temperature (15-25°C).

## PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

## PROCEDURE

Wavelength: 505 nm (490 - 550)  
Lightpath: 1 cm  
Temperature: 37°C / 15-25°C  
Reading: against distilled water

pipette:	blank	sample	standard
Reagent (A+B)	1000 µl	1000 µl	1000 µl
water	10 µl		
sample		10 µl	
standard			10 µl

Mix, incubate for 5 minutes at 37°C or 10 minutes at room temperature (15-25°C).

Read against blank the absorbance of the samples (Ax) and the standard (As).  
The color is stable at least 30 minutes.

## RESULTS CALCULATION

**Lactate (mg/dl)** =  $A_x / A_s \times 10$  (Standard conc.)

Conversion Factor: mg/dl  $\times 0.111$  = mmol/l

## EXPECTED VALUES

**4.5 – 19.8 mg/dl (0.5 – 2.2 mmol/l)**

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

## PERFORMANCE

**Sensitivity:** the sensitivity of the method is: 0.39 mg/dl.

**Linearity:** the linearity is up to 150 mg/dl. For higher values, dilute the sample 1:2 and multiply the result by 2.

## Precision intra-assay:

	Level 1	Level 2
Mean (mg/dl)	11.1	21.8
DS	0.24	0.25
CV %	2.14	1.16

## Precision inter-assay:

	Level 1	Level 2
Mean (mg/dl)	11.4	22.1
DS	0.36	0.54
CV %	3.12	2.47

**Interferences:** Intravenous injection of epinephrine, glucose, bicarbonate, or other infusions that modify the acid-base balance, cause an elevation in lactate.

**Correlation against a reference method:**  $Y = 0.9979x - 1.2518$   $r = 0.998$ .

## NOTE

Calibration with the aqueous standard may cause a systematic error in automatic procedures. In these cases, it is recommended to use a serum Calibrator.

## REFERENCES

- Gau N. Lactic acid. Kaplan A et al. Clin Chem The C.V. Mosby Co. St. Louis. Toronto. Princeton 1984; 1040-1042 e 418.
- Young D. S., Effects of drugs on Clinical Lab. Tests, . 4th ed. AACC Press (1995).
- Young D. S., Effects of disease on Clinical Lab. Tests, 4th ed. AACC (2001).
- Burtis A. et al. Tietz Textbook of Clinical Chemistry, 3th ed. AACC (1999).
- Tietz N. W. et al. Clinical Guide to Laboratory Tests, 3th ed. AACC (1995).