# LIPASI

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# LIPASE

Colorimetric Method  
Liquid Reagent ready to use

REF. 5722 4x10+1x10 ml

DNV CERTIFIED COMPANY

UNI EN ISO 9001:2008  
UNI EN ISO 13485:2012

## INTENDED USE

Quantitative determination of Lipase in serum and plasma.

## PRINCIPLE

The colorimetric substrate, 1,2-O-Dilauryl-rac-glycero-3-glutaric acid-(6'-methyl-resorufin)-ester, is cleaved by pancreatic lipase and the resulting dicarboxylic acid ester is hydrolysed under the alkaline test conditions to yield the chromophore methylresorufin.

The formed methylresorufin, photometrically measured, is proportional to the amount of lipase in the sample.

## SAMPLE

Serum, plasma with heparin.

Lipase activity is stable 7 days in samples stored at 2-8 °C.

## KIT COMPONENTS

Reagent (A) Volume = 10 ml	Good's Buffer pH 8.0, colipase $\geq 2$ mg/l, desoxycolate $\geq 1.0$ mM, taurodesoxycolate $\geq 1.0$ mM, calcium ions $\geq 1$ mM, detergent and preservative.
Reagent (B) Volume = 10 ml	Tartrate buffer pH 4.0, lipase substrate $\geq 0.1$ mM, stabilizer and preservative.

The reagents are stable until the expiration date indicated on the label if stored at 2-8°C. Do not use over expiry date.

Once opened reagents are stable for 2 months at 2-8°C if contamination is avoided. Keep bottles closed when not in use.

The Reagent (B) is an orange-colored micro-emulsion. A slight apparent precipitation could occur, showing a light red deposit on the bottom of vial. It is a normal behaviour and it is recommended to resuspend solution before analysis, with a mild shaking.

## REAGENT PREPARATION

Liquid Reagents, ready to use.

Reagent (B): Mix gently before use.

## PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

## PROCEDURE

Wavelength: 580 nm (570-590 nm)

Lightpath: 1 cm

Temperature: 37°C

Reading: against distilled water

pipette:	blank	sample	calibrator
Reagent (A)	1000 $\mu$ l	1000 $\mu$ l	1000 $\mu$ l
water	20 $\mu$ l		
sample		20 $\mu$ l	
calibrator			20 $\mu$ l

Mix carefully, incubate at 37°C for 5 minutes.

Reagent (B)	250 $\mu$ l	250 $\mu$ l	250 $\mu$ l
-------------	-------------	-------------	-------------

Mix, execute a first reading of absorbance after 2 minutes, incubating at 37°C. Perform other 2 readings at 60 seconds intervals. Calculate the  $\Delta A/\text{min}$ .

## RESULTS CALCULATION

$$\Delta A/\text{min} = \Delta A/\text{min}_{(\text{sample or calibrator})} - \Delta A/\text{min}_{(\text{blank})}$$

$$\text{Lipase (U/L)} = [\Delta A/\text{min}_{(\text{sample})} / \Delta A/\text{min}_{(\text{calibrator})}] \times \text{Calibrator Value}$$

## EXPECTED VALUES

$\leq 60.0$  U/L

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

## PERFORMANCE

**Sensitivity:** the sensitivity of the method is 1 U/l.

**Linearity:** the method is linear up to 300 U/l. For higher values, dilute the sample 1:2 and multiply the result by 2.

**Precision intra-assay (n=10):**

	Level 1	Level 2
Mean (U/l)	60.6	90.4
DS	0.54	0.70
CV %	0.89	0.77

**Precision inter-assay (n=20):**

	Level 1	Level 2
Mean(U/l)	59.9	90.3
DS	1.76	1.80
CV %	2.94	1.99

**Interferences:** bilirubin does not interfere up to 50 mg/dl. Hemoglobin does not interfere up to 400 mg/dl. Ascorbic acid does not interfere up to 50 mg/dl. Lipids do not interfere up to 1000 mg/dl.

**Correlation against a reference method:**  $Y = 0.93x + 0.50$   $r = 0.99$ .

## REFERENCES

1. Tietz N. and Shuey DF., Clin. Chem. 1993, 39, 746-756.
2. Tietz Textbook of Clinical Chemistry, Fourth Edition, Burtis-Ashwood-Brunts (2006) 619-621.

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) - web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

572207

Ed. 2018/06 rev. 00