# POTASSIO

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Ed. 2018/03 rev. 02

# POTASSIUM

Method: Enzymatic  
Liquid Reagents, ready to use

REF. 0236 2x50 ml + Cal 2x2 ml  
REF. 0236/2 1x50 ml + Cal 2x1 ml

DNV CERTIFIED COMPANY  
UNI EN ISO 9001:2008  
UNI EN ISO 13485:2012

## INTENDED USE

Quantitative determination of potassium in serum.  
Measurements obtained by this assay are used to monitor electrolyte balance in the diagnosis and treatment of diseases conditions characterized by low or high blood potassium levels.

## PRINCIPLE

Potassium is determined spectrophotometrically through a kinetic coupling assay system using potassium dependent pyruvate kinase. Pyruvate generated is converted to lactate accompanying conversion of NADH analog to NAD analog. The corresponding decrease of optical density at 380 nm is proportional to the potassium concentration in the serum.

## SAMPLE

Non-haemolysed serum  
It is advisable to collect the samples so you can test them as soon as possible and within five days from samples collection.

## KIT COMPONENTS

Reagent (A) Potassium Volume = 2x40 ml	LDH (< 50 KU/L), Analog (< 10 mmol/L), NADH, Substrate, Azide(0.05 %), Stabilizers.
Reagent (B) Potassium Volume = 2x10 ml	Pyruvate kinase (< 50 KU/L), Azide (0.05 %), Stabilizers.
Low Calibrator Volume = 1-2 ml	Potassium (batch-specific value indicated on the bottle)
High Calibrator Volume = 1-2 ml	Potassium (batch-specific value indicated on the bottle)

Reagents and calibrators are stable until the expiration date indicated on the label if stored at 2-8°C. Do not use over the expiry date.  
Once opened reagents are stable for 2 months if contamination is avoided.  
Keep bottles closed when not in use.

## REAGENT PREPARATION


Liquid Reagents, ready to use.

## PRECAUTIONS AND WARNINGS

Use the normal precautions required in the laboratory.  
Do not swallow. Avoid contact with eyes and skin.  
These reagents contain lithium azide. The azide may react with lead or copper plumbing to form potentially explosive compounds. Upon disposal of these reagents, flush drains with water. Dispose of waste in accordance with local regulations.  
The Material Safety Data Sheet is available on request.

## PROCEDURE

Diagram of the test to automated chemistry analyzers:

## CALIBRATION

For the calibration you must use the standards high and low for potassium supplied. The concentration of the potassium in the sample is determined by a linear calibration curve using the standards high and low.

In the case of analyzers that require a zero calibration, it is possible to use a saline solution for this purpose and use standard high and low potassium provided as calibrators respectively 2 and 3.  
It is advisable to calibrate weekly.

## REFERENCE RANGE

The normal level of potassium in human serum is in a range between:

**3.5 – 5.1 mmol/L (13.7 - 19.9 mg/dl)**

Each laboratory should establish appropriate reference intervals related to its population.

## QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions.

## PERFORMANCE

**Linearity:** the method is linear for concentrations of potassium between 2.0 and 8.0 mmol/L.

**Limit of detection:** the detection limit of the enzymatic test for potassium is equal to 0.87 mM K<sup>+</sup>.

**Precision intra-assay:**

	Level 1	Level 2
Mean (mM)	4.62	6.96
CV %	1.12	1.20

**Precision inter-assay:**

	Level 1	Level 2
Mean (mM)	4.62	6.96
CV %	1.77	1.77

**Interferences:** the test is not affected by the interference of the following substances, if present at the indicated concentrations: Na<sup>+</sup> to 150 mM, NH<sub>4</sub><sup>+</sup> to 0.5 mM, Ca<sup>2+</sup> to 7.5 mM, Pi to 2.0 mM, ascorbic acid to 10.0 mM, Zn<sup>2+</sup> to 0.5 mM, Fe<sup>3+</sup> to 0.5 mM, Cu<sup>2+</sup> to 0.5 mM, haemoglobin to 500 mg/dl, triglycerides to 1000 mg/dl, conjugated bilirubin to 20 mg/dl and unconjugated bilirubin to 15 mg/dl.

**Correlation against a reference method:** Y = 1.07x - 0.30 r = 0.98

## REFERENCES

- Wu, A.H.B., ed. Tietz Clinical guide to Laboratory Tests, 4<sup>th</sup> ed., p. 880. W.B. Saunders Company, St. Louis (2006).
- Bergmeyer, H.U., Gawehn, K., and Grassl, M. (1974) in *Methods of Enzymatic Analysis*. Second Edition, Volume I, 509-510, Academic Press, Inc., New York.
- M.N. Berry, R.D. Mazzachi, M. Pejakovic, and M.J. Peake Enzymatic Determination of Potassium in Serum. CLIN.CHEM. 35/5, 817-820 (1989).