PROTEINE TOTALI

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TOTAL PROTEINS

Colorimetric Method - Biuret
Liquid Reagent ready to use

REF. 0046/2 2x 100 ml
REF. 0046 4x 100 ml
REF. 0044 2x 250 ml
REF. 0042 4x 250 ml
REF. 0043 1x1000 ml

DNV CERTIFIED COMPANY
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

INTENDED USE

Quantitative determination of total proteins in serum and plasma.

PRINCIPLE

In an alkaline medium, proteins form a blue-violet complex with Cu (II) ions. The color intensity is proportional to the concentration of total proteins in the sample.

SAMPLE

Serum, plasma with heparin. Do not use hemolyzed samples.
Proteins in the sample are stable one week at room temperature (15-25 °C) and one month at 2-8°C.

KIT COMPONENTS

| | | |
|--|--|--|
| Reagent (A) PT Volume = 100/250/1000 ml | Sodium hydroxide Copper sulphate Potassium iodide Potassium sodium tartrate | 200 mmol/l 8 mmol/l 33 mmol/l 10 mmol/l |
|--|--|--|

Optional: Standard for Total Proteins REF. 8801 o 8801/6

The Standard is not included in the kits.

| | | |
|----------------------------------|--|---------------------|
| Standard (B) PT Volume = 5 ml | Total Proteins (Albumin) Sodium azide | 6 g/dl 14 mmol/l |
|----------------------------------|--|---------------------|

The reagents are stable until the expiration date if stored at temperatures indicated on the label and protected from light. Once opened, the reagent are stable 2 months if contamination is avoided.

Keep bottles closed when not in use.

REAGENT PREPARATION

Liquid Reagent ready to use.

PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

Warning

- H319** Causes serious eye irritation.
H315 Causes skin irritation.
H411 Toxic to aquatic life with long lasting effects..
- EUH208** Contains: POTASSIUM IODIDE
May produce an allergic reaction.
- P264** Wash with water thoroughly after handling.
P273 Avoid release to the environment.
P280 Wear protective gloves / protective clothing/ eye protection/ face protection.
- P302+P352** IF ON SKIN: Wash with plenty of soap and water.
P332+P313 If skin irritation occurs: Get medical advice / attention.

Reagents may contain some non-reactive and preservative components. Use the normal precautions required in the laboratory.
Dispose of waste according to local laws.

PROCEDURE

Wavelength: 546 nm (530 – 550)
Lightpath: 1 cm
Temperature: 25, 30, 37°C
Reading: against blank reagent
Method: Increasing End Point
Sample/Reagent: 1/100

| pipette: | blank | Sample | standard |
|-------------|---------|---------|----------|
| Reagent (A) | 1000 µl | 1000 µl | 1000 µl |
| water | 10 µl | | |
| sample | | 10 µl | |
| standard | | | 10 µl |

Mix, incubate at 25, 30, 37°C for 10 minutes and read against blank reagent the absorbance of the sample (Ax) and the standard (As).

Reaction volumes can be proportionally varied.
This method describes the manual procedure to use the kit.
For automated procedure, ask for specific applications.

RESULTS CALCULATION

Total Proteins g/dl = $A_x/A_s \times 6$ (standard Value)

EXPECTED VALUES

Serum/plasma: 6.6 – 8.3 g/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL

You must perform the controls at each kit's use and verify that the values obtained are within the reference range reported in the operating instructions. For this purpose we recommend the use of control sera: PRECISENORM (REF.6000) and PRECISEPATH (REF.6001).

PERFORMANCE

Sensitivity: the sensitivity of the method is: 0.1 g/dl.

Linearity: the method is linear up to 13 g/dl. For higher values, dilute with saline 1:2 the sample and multiply the result by 2.

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| | Level 1 | Level 2 | Level 3 |
|-------------|---------|---------|---------|
| Mean (g/dl) | 2.58 | 5.08 | 9.75 |
| DS | 0.043 | 0.044 | 0.110 |
| CV % | 1.65 | 0.87 | 1.13 |

Precisione tra le serie:

| | Level 1 | Level 2 | Level 3 |
|-------------|---------|---------|---------|
| Mean (g/dl) | 2.65 | 5.36 | 10.12 |
| DS | 0.028 | 0.052 | 0.156 |
| CV % | 1.04 | 0.97 | 1.54 |

Interferences: bilirubin does not interfere up to 30 mg/dl. A triglycerides concentration exceed 300 mg/dl increases the reading.

Hemoglobin in concentration > 0.3 mg/dl increases the reading.

Correlation against a reference method: $Y = 1.0068x - 0.0670$ $r = 0.9974$

REFERENCES

1. Vassault, A. et al. Ann. Biol. Clin., 44,686 (1986).
2. Young, D.S., et al., Clin. Chem. 21:1D (1975).
3. Kaplan LA, Pesce AJ, Clinical Chemistry, Mosby Ed. 1989
4. Kingsley G.R., Biol. Chem. 131, 197-200 (1939)

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18,300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

004607

Ed. 2014/12 rev. 03