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## BIBLIOGRAFIA

1. NCCLS: One - Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; approved Guideline. NCCLS document H47 - A, NCCLS, PA, (1996).
2. NCCLS Collection, Transport and Processing of blood Specimens for Coagulation Testing and Performance of Coagulation Assays 2<sup>nd</sup> Ed. Approvate guideline. NCCLS Document H21 - A3. Wayne, PA (1998).

## INTENDED USE

This Calibrator is a human plasma especially designed to build standard curves. It is a valid instrument which guaranties an optimum calibration on automatic analyzers.

## KIT COMPONENTS

Calibrator ( 1 ml)	Human Plasma collected with < 0.4 % sodium citrate, Hepes 0.6 %, NaOH 0.08%, Lithium Chloride 0.04 %, Polybrene 0.001%.
-----------------------	---

Values of each test is carried in the list present in the package.

## PREPARATION

Reconstitute Lyophilized Calibrator with 1 ml of distilled water. Mix gently and leave stand at room temperature for 15 minutes (Do not invert vial or mix vigorously).

## STABILITY AND STORAGE

Store vials tightly closed at 2-8°C. Stable until the expiration date stated on the label.

After proper reconstitution, the control is stable 8 hours at 2-8°C.

**Bring the Calibrator at room temperature before use.**

## PRECAUTIONS AND WARNINGS

Whenever infectious agents, chemical reagents, reagents of human/animal origin, blood or other biological liquids are manipulated, it is advisable to follow the most common recommendations and take all the necessary hygienic precautions as the monouse gloves.

Each unit of source material used in the preparation of Calibrator has been tested by an approved method and found non-reactive for HbsAg, and negative for antibodies to HCV and HIV 1/2.

However, no Known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS or other infectious diseases. This product, like all materials of human origin, should be handled as potentially infectious material.

Dispose of waste according to local laws.

## RESULTS

Results depend on reagents, instruments and individual laboratory protocols.

## LIMITATIONS

Calibrator is subject to the limitations of the test system.

Variables such as temperature, reagent stability, instrument performance and individual technique can influence final results.

Always follow instrument and reagent manufacturers guideline.

## NOTES

Sometimes the Calibrator have a slight turbidity. It is right to eliminate the product only in case of excessive turbidity, bacterial contamination or colour change.

## REFERENCES

1. NCCLS: One - Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; approved Guideline. NCCLS document H47 - A, NCCLS, PA, (1996).
2. NCCLS Collection, Transport and Processing of blood Specimens for Coagulation Testing and Performance of Coagulation Assays 2<sup>nd</sup> Ed. Approvate guideline. NCCLS Document H21 - A3. Wayne, PA (1998).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) - web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

101207  
Ed. 2012/10 rev. 00