CALBRATORE FIBRINOGENO

REF. 1006 5x1 ml

AZIENDA CERTIFICATA DNV
UNI EN ISO 9001:2008
UN EN ISO 13485:2012

USO PREVISTO

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BIBLIOGRAFIA

NCCLS: Procedure for determining fibrinogen in plasma. Approved Guideline. NCCLS document H30 – A, Villanova, PA, (1994).

INTENDED USE

The Fibrinogen Calibrator is a human lyophilized plasma intended for use as Calibrator in quantitative assay of Fibrinogen and ensure optimum calibration for use in automatic instruments.

KIT COMPONENTS

Calibrator (1 ml)	Human Plasma collected with < 0.1 % sodium citrate, Hepes 0.4 %, NaOH 0.08%, Lithium Chloride 0.04 %, Polybrene 0.001%.
-------------------	---

Values is carried on the label.

PREPARATION

Reconstitute Lyophilized Calibrator with 1 ml of distilled water.

Mix gently and leave stand at room temperature for 15 minutes (Do not invert vial or mix vigorously).

STABILITY AND STORAGE

Store vials tightly closed at 2-8°C. Stable until the expiration date stated on the label.

After proper reconstitution, the calibrator is stable 8 hours at 2-8°C.

Bring the Calibrator at room temperature before use.

PRECAUTIONS AND WARNINGS

Whenever infectious agents, chemical reagents, reagents of human/animal origin, blood or other biological liquids are manipulated, it is advisable to follow the most common recommendations and take all the necessary hygienic precautions as the monouse gloves.

Each unit of source material used in the preparation of Control has been tested by an approved method and found non-reactive for HbsAg, and negative for antibodies to HCV and HIV 1/2.

However, no Known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS or other infectious diseases.

This product, like all materials of human origin, should be handled as potentially infectious biological material.

Dispose of waste according to local laws.

TEST PROCEDURE

After reconstitution, prepare a minimum of 5 different dilutions with Imidazole buffer: 1:5, 1:10, 1:15, 1:20, 1:40.

To use the dilutions as samples.

Use all five of the calibrator points to construct a log-log curve that plots fibrinogen concentration (mg/dl) vs. clotting time.

Derive the concentration of fibrinogen in the sample reporting on the calibration curve its clotting time.

RESULTS

For each lot value will given.

Results depends in great amount on calibrator, reagent, instruments and individual laboratory protocols.

Each laboratory should make own evaluation of product on the grounds of reagents, instruments and protocols to establish own reference range.

NOTE

Sometimes the Calibrators have a slight turbidity.

It is right to eliminate the product only in case of excessive turbidity, bacterial contamination or colour change.

LIMITATIONS

Calibrators are subject to the limitations of the test system.

Variables such as temperature, reagent stability, instrument performance and individual technique can influence final results.

Always follow instrument and reagent manufacturers guideline.

REFERENCES

NCCLS: Procedure for determining fibrinogen in plasma. Approved Guideline. NCCLS document H30 – A, Villanova, PA, (1994).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

100607

Ed. 2014/12 rev. 01