# CONTROL PLASMA P

REF. 1008

5x1 ml

AZIENDA CERTIFICATA DNV

UNI EN ISO 9001:2008

UN EN ISO 13485:2012

## USO PREVISTO

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## BIBLIOGRAFIA

1. NCCLS: One - Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; approved Guideline. NCCLS document H47 - A, NCCLS, PA, (1996).
2. NCCLS Collection, Transport and Processing of blood Specimens for Coagulation Testing and Performance of Coagulation Assays 2<sup>nd</sup> Ed. Approvate guideline. NCCLS Document H21 - A3. Wayne, PA (1998).

## INTENDED USE

Coagulation Controls are intended for use as control to monitor the performance of routine coagulation assays.

## KIT COMPONENTS

Control (5x1 ml)	Human Plasma collected with < 0.4 % sodium citrate Hepes 0.4 %, NaOH 0.3%, Lithium Chloride 0.08%, Barium sulphate 0.1%, Aprotinin 0.01%, Polybrene 0.001%.
------------------	---

Values of each test is carried in the list present in the package.

## PREPARATION

Reconstitute Lyophilized Control with 1 ml of distilled water. Mix gently and leave stand at room temperature for 15 minutes (Do not invert vial or mix vigorously).

## STABILITY AND STORAGE

Store vials tightly closed at 2-8°C. Stable until the expiration date stated on the label.

After proper reconstitution, the control is stable 8 hours at 2-8°C.

Bring the Control at room temperature before use.

## PRECAUTIONS AND WARNINGS

Whenever infectious agents, chemical reagents, reagents of human/animal origin, blood or other biological liquids are manipulated, it is advisable to follow the most common recommendations and take all the necessary hygienic precautions as the monouse gloves.

Each unit of source material used in the preparation of Control has been tested by an approved method and found non-reactive for HbsAg, and negative for antibodies to HCV and HIV 1/2.

However, no Known test method can offer complete assurance that products derived from human blood will not transmit hepatitis, AIDS or other infectious diseases.

This product, like all materials of human origin, should be handled as potentially infectious biological material.

Dispose of waste according to local laws.

## TEST PROCEDURE

Use Control as samples.

See the procedure of used reagent.

## RESULTS

Results depend on reagents, instruments and individual laboratory protocols.

Values are expected to be moderately abnormal.

## PERFORMANCE

Control Plasma P have been adjusted to have prolonged prothrombin (PT) and partial thromboplastin (APTT) times.

Actual results depend on lot number, type of reagent and instrument.

## LIMITATIONS

Control is subject to the limitations of the test system.

Variables such as temperature, reagent stability, instrument performance and individual technique can influence final results.

Always follow instrument and reagent manufacturers guideline.

## NOTES

Sometimes the Control have a slight turbidity. It is right to eliminate the product only in case of excessive turbidity, bacterial contamination or colour change.

## REFERENCES

1. NCCLS: One - Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; approved Guideline. NCCLS document H47 - A, NCCLS, PA, (1996).
2. NCCLS Collection, Transport and Processing of blood Specimens for Coagulation Testing and Performance of Coagulation Assays 2<sup>nd</sup> Ed. Approvate guideline. NCCLS Document H21 - A3. Wayne, PA (1998).

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giesseiagnostics.com](mailto:info@giesseiagnostics.com) - web site: [www.giesseiagnostics.com](http://www.giesseiagnostics.com)

100807

Ed. 2014/12 rev. 01