
# VDRL

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REF. 7800 - 250 tests

AZIENDA CERTIFICATA DNV  
UNI EN ISO 9001:2008  
UN EN ISO 13485:2012

## USO PREVISTO

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## BIBLIOGRAFIA

1. George P. Schmid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra Larsen et al. Clinical Microbiology Reviews 1995; 8 (1) : 1-21.
3. Sandra Larsen et al. A manual test for Syphilis American Public Health Association 1990: 1-192.
4. Young DS. Effects of drugs on Clinical Laboratory Tests, 4th ed. AACC Press 1995.

Giesse Diagnostics srl

V. Enrico Fermi,3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111


e-mail: [info@giesseiagnostics.com](mailto:info@giesseiagnostics.com) - web site: [www.giesseiagnostics.com](http://www.giesseiagnostics.com)

780007

Ed. 2014/12 rev. 01

# VDRL

Slide agglutination

REF. 7800 - 250 tests

DNV CERTIFIED COMPANY

UNI EN ISO 9001:2008  
UNI EN ISO 13485:2012

## INTENDED USE

Qualitative determination of plasma reagins.

## PRINCIPLE

The VDRL test is a non-treponemal slide agglutination test for the qualitative and semi-quantitative detection of plasma reagins.

The antigen suspension, a lipid complex, is agglutinated when mixed with samples containing reagins of patient affected by syphilis.

## SAMPLE

Fresh Serum, plasma or cerebrospinal fluid. Stable 7 days at 2-8°C or 3 months at -20°C.

Samples with presence of fibrin should be centrifuged before testing.

Do not use highly hemolized or lipemic samples.

## KIT COMPONENTS

Reagent (A) VDRL antigen Volume = 5.0 ml	Solution containing cardiolipin 0.3 g/L, lecithin 2.1 g/L and cholesterol 9 g/L in phosphate buffer 1.5 mmol/l. Preservative. pH 7.0
Control (+) VDRL Volume = 1.0 ml	Artificial serum with reagin titer $\geq 1/8$
Control (-) VDRL Volume = 1.0 ml	Animal serum. Preservative

The Reagents are stable until the expiration date printed on the label, when stored tightly closed at 2-8°C. Once opened, the reagents are stable one month at 2-8°C if contamination is avoided. Do not freeze.

Keep bottles closed when not in use.

## PRECAUTIONS AND WARNINGS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.

Use the normal precautions required in the laboratory.

Dispose of waste according to local laws.

## REAGENT PREPARATION

The Reagents are ready to use.

## ADDITIONAL EQUIPMENT

Mechanical rotator with adjustable speed at 180 r.p.m.

Glass slides.

Light microscope (100x).

## PROCEDURE

### Qualitative method:

Allow the reagents and samples to reach room temperature. The sensitivity of the test may be reduced at low temperatures.

Place 50  $\mu$ l of the sample and one drop of each Positive and Negative Controls into separate circles on the slide test.

Swirl the VDRL suspension gently before using and add 20  $\mu$ l of this reagent onto each sample.

Place the slide on a mechanical rotator at 180 rpm for 4 minutes.

False positive results could appear if the test is read later than 4 minutes.

### Semi-quantitative method:

Make serial two fold dilutions of the sample in saline solution.

Proceed for each dilution as in the qualitative method.

## READING AND INTERPRETATION

Examine the presence or absence of agglutination immediately after rotation using the light microscope (100x).

Agglutination	Report
Medium or large clumps	Reactive
Small clumps	Weakly reactive
No clumping or very slight "roughness"	Non Reactive

The titer, in the semi-quantitative method, is defined as the highest dilution showing a positive result.

## QUALITY CONTROL

Positive and Negative Controls are recommended to monitor the performance of the reagent and to have a better results interpretation.

## PERFORMANCE

**Prozone effect:** No prozone effect up to titers  $\geq 1/128$

**Diagnostic sensitivity:** 100 %

**Diagnostic specificity:** 100 %

**Interferences:** bilirubin does not interfere up to 20 mg/dl. Lipids and hemoglobin do not interfere up to 10 g/l. Rheumatoid factors ( 300 U/ml ) interfere. Other substances may interfere.

## METHOD LIMITATIONS

VDRL test is non-specific for syphilis. All Reactive samples should be retested with treponemal methods such as TPHA to confirm the results.

A Non Reactive result by itself does not exclude a diagnosis of syphilis.

False positive results have been reported in diseases such as infectious mononucleosis, viral pneumonia, toxoplasmosis, pregnancy and autoimmune diseases.

## REFERENCES

1. George P. Schimid. Current Opinion in Infectious Diseases 1994; 7: 34-40.
2. Sandra Larsen et al. Clinical Microbiology Reviews 1995; 8 (1) : 1-21.
3. Sandra Larsen et al. A manual test for Syphilis American Public Health Association 1990: 1-192.
4. Young DS. Effects of drugs on Clinical Laboratory Tests, 4th ed. AACC Press 1995.

## Giesse Diagnostics srl

V. Enrico Fermi,3 - Z.I. V. Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: [info@giessediagnostics.com](mailto:info@giessediagnostics.com) - web site: [www.giessediagnostics.com](http://www.giessediagnostics.com)

780007

Ed. 2014/12 rev. 01