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INTENDED USE

The systemic solution is a solution that can be used in all automated instruments to allow a better sliding of the liquid in the hydraulic circuit of the analyzer, and thus prevent the adhesion of small drops of liquid on the walls of the tubes, on the needle, on the reading cuvette.

The new formulation avoids the pollution between reagents and does not interfere with the reagents used in the tests.

PRODUCTS CHARACTERISTICS

Solution (A): Surfactant with antibacterial.

Solution (B): HCl Solution.

PREPARATION

Use each solution in the ratio: 1:1001.

Pour 10 ml of Solution (A) and 10 ml of Solution (B) in the container of the washing water of the analyzer and add 10 liters of distilled water.

You should prepare a weekly wash solution, eliminating the old residual solution and thoroughly rinse the container with distilled water.

STORAGE AND STABILITY

The product is stable until the expiration date indicated on the label when stored at 2-25°C, protected from light and in the absence of contamination.

Keep bottles closed when not in use.

PRECAUTIONS AND WARNINGS

EC Regulation 1272/2008 (CLP):

Soluzione (B):

Warning

- H319** Causes severe eye irritation.
H315 Causes skin irritation.
H335 May cause respiratory irritation.
P264 Wash thoroughly with water after use.
P280 Wear protective gloves/ protective clothing/ eye protection/ face
P312 If feel unwell, contact a POISON CENTER or physician.
P304+P340 If Inhaled: remove victim to fresh air and keep at rest in a position comfortable for breathing.
P403+P233 keep container tightly closed and in a well-ventilated area.

Use in accordance with good laboratory practise, avoiding contact with skin, eyes and swallow.

WASTE DISPOSAL

Dispose of product according to local laws.

Giesse Diagnostics srl

V. Enrico Fermi, 3 - Z.I. V.Tiburtina Km 18.300 - 00012 Guidonia Montecelio (RM) - Italia

Tel. +39 0774 051100 - Fax +39 0774 051111

e-mail: info@giessediagnostics.com – web site: www.giessediagnostics.com

Ed. 2016/01 rev. 05